



POLYPRO PREFABS

Pre-stressed, semi rigid polypropylene prefabricated devices that incorporate many features of custom orthoses at a fraction of the price.

Polypro Prefabs are ideal as an economical first step in the treatment of a variety of common foot, ankle and lower leg pathologies that are secondary to mild to moderate subtalar joint pronation, excessive mid-tarsal joint motion and/or mild forefoot deformities.

Polypro Prefabs are also an excellent choice as an inexpensive "second pair" for those patients with limited insurance coverage.



Polypro Prefabs were designed by ProLab Orthotics/USA, one of America's most respected prescription foot orthoses labs, and incorporate a number of advanced functional control features including:

- 2mm medial heel skive to enhance rear foot control and decelerate subtalar joint pronation,
- 3 degree intrinsic forefoot valgus correction to load the lateral column and limit mid-tarsal joint motion,
- Traditional "functional" shell shape parameters,
- 14mm heel cup, and
- A textured surface that eliminates the need for a top cover.

Polypro Prefabs are sized by shoe size and are available in men's and women's sizes 6 through 13. Polypro Prefabs are available through foot healthcare professionals' offices across Canada.



PROLAB
ORTHOTICS/USA

Paris orthotics