

Paris

## WALKWELL 328 UNIVERSAL POST OPERATIVE SHOE



**WalkWell's new 328 Universal Post Operative Shoe** is a welcome alternative to the many poorly constructed, ill-fitting post-op shoes currently on the market.

The **WalkWell 328** Features

- Rigid rocker sole that effectively offloads metatarsal head pressure by eliminating metatarsal-phalangeal joint dorsiflexion,
- Bacteria static EVA insole that comfortably attenuates shock,
- Full reinforced heel counter with a unique adjustable closure system designed to accommodate heels of almost any size or shape,
- Soft and breathable micro-fabric upper with three Velcro closures to gently accommodate edema, post-op dressings or accommodative orthoses.
- Aesthetically pleasing all black design that promotes patient compliance.

### Clinical Indications

- Osteotomies with and without fixation
- Post nail surgery
- Metatarsal fractures
- Metatarsal phalangeal joint injuries
- Metatarsal phalangeal synovitis and tendonitis
- Hallux limitus
- Hallux Valgus
- Neuropathic ulceration
- Persistent edema

The **WalkWell 328** rigid rocker sole design is the best in its class. A study published in the Journal of the American Podiatric Medical Association indicates that WalkWells can reduce peak forefoot pressures up to 20% better than conventional post-op shoes.

**WalkWell 328's** are universal (fit left or right) and are available in seven sizes from Women's Small (5) to Men's XX-Large (16.5). 328's are sold in single units as well as a 14 shoe (mixed size) starter package.

For more information or to request an order form please call:

**Paris** orthotics

Paris Orthotics Ltd.  
3630 East 1<sup>st</sup> Avenue  
Vancouver BC V5M 1C3  
Toll Free: 1.800.848.0838

