



PRE FABRICATED ORTHOTICS

Paris Orthotics provides over the counter orthotics for the whole family through our Adult Prefab Orthotics and Kiddythotics.



SHIPPING DETAILS

SHIP VIA	
<input type="checkbox"/>	Include this order with my next orthotic shipment at no extra charge
<input type="checkbox"/>	Ship ASAP via UPS - bulk shipping charge will be added to your order
SHIP TO	
Account Name _____	
Clinician _____	
Phone _____	

KIDDYTHOTICS

An easy and effective over the counter option for the paediatric patient. Each Kiddythotic Prefab is colour coded by size making orthotics fun for children and an easy way to distinguish one from another.

Each device is made with a rigid polypropylene shell which is durable, resistant to compression and provides consistent control. The medial skive supports the medial column and attenuates pronation, while the rearfoot post and medial flange provide rearfoot and midfoot stability. Lastly, the deep heel cup improves calcaneal stabilization. The Kiddythotic Prefab is a great over the counter and economical solution for the rapidly growing child and youth population.

- POLYPROPYLENE SHELL
- MEDIAL FLANGE
- DEEP HEEL CUP
- 4/4 REARFOOT POST
- 4MM HEEL SKIVE



PRODUCT CODE	ESTIMATED SHOE SIZE	ORDER QUANTITY	PRICE PER PAIR
A Purple			
B Yellow			
C Orange			
D Green			
E Red			
F Blue			
G Grey			
Sizing Set			

6+ pair orders receive a complimentary sizing set.



PRE FABRICATED ORTHOTICS

Paris Orthotics provides over the counter orthotics for the whole family through our Adult Prefab Orthotics and Kiddythotics.



ADULT PREFAB

An Adult Prefab is a great solution for patients who need immediate support. Adult Prefabs are available in two options; unposted (natural) and posted (black). Both over the counter options are well suited for dress casual footwear where the insole is not removable.

Each is made with a semi rigid polypropylene shell which is durable, resists compression and ensures control. The medial skive supports the medial column to attenuate pronation. The valgus post stabilizes the midtarsal joint to reduce forefoot motion addressing the valgus forefoot deformity found in most pathological feet. The difference between the two inserts is the rearfoot post option; the posted device includes a rearfoot post which provides increased rearfoot and midfoot control. The Adult Prefabs are an excellent choice for patients who would benefit from any of the features above.

SHIPPING DETAILS

SHIP VIA

Include this order with my next orthotic shipment at no extra charge

Ship ASAP via UPS - bulk shipping charge will be added to your order

SHIP TO

Account Name _____

Clinician _____

Phone _____

- POLYPROPYLENE SHELL
- 14MM HEEL CUP DEPTH
- 2MM HEEL SKIVE
- 3 DEGREES FOREFOOT VALGUS CORRECTION
- REARFOOT POST (POSTED ONLY)



ORTHOTIC SIZE		POSTED QUANTITY	PRICE PER PAIR	UNPOSTED QUANTITY	PRICE PER PAIR
W-6					
W-7					
W-8	M-6				
W-9	M-7				
W-10	M-8				
W-11	M-9				
W-12	M-10				
W-13	M-11				
	M-12				
	M-13				

10+ pair orders receive a 10% discount.