

# DIABETIC & RHEUMATOID ARTHRITIS PREFAB ORTHOSES

## A COST EFFECTIVE OPTION FOR AT RISK FEET

This line of prefabricated insoles offers unique treatment solutions for a variety of control and protection needs often present in the Diabetic and Rheumatoid Arthritic populations. These prefab orthoses are available in two different designs; SOFT and GENTLE and can be an economical solution to provide interim protection and support while custom orthoses are being fabricated.

### SOFT

#### DRA | 1

An off the shelf insole for those that require maximum cushioning and protection.



#### FEATURES

- Smooth Poron Top Cover
- 15 Durometer Polyurethane Insole
- Available in Shoe Sizes from Women's 5 to Men's 15

#### APPLICATIONS

- At risk feet that require exceptional cushioning and outstanding pressure relief
- A first response treatment prior to dispensing custom foot orthoses

### GENTLE

#### DRA | 2

An off the shelf insole for those that require mild support, cushioning and protection.



#### FEATURES

- P-Cell top cover
- 40 durometer EVA insole
- Available in shoe sizes from Women's 5 to Men's 15

#### APPLICATIONS

- To use when mild control and shear force diffusion is needed
- A first response treatment prior to dispensing custom foot orthoses

**DIABETIC & RHEUMATOID ARTHRITIS PREFAB ORTHOSES**
**SHIP TO**

ACCOUNT NAME \_\_\_\_\_

ACCOUNT # \_\_\_\_\_

CLINICIAN \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY \_\_\_\_\_

PROVINCE \_\_\_\_\_ POSTAL CODE \_\_\_\_\_

PHONE \_\_\_\_\_ FAX \_\_\_\_\_

EMAIL \_\_\_\_\_

**ORDER SUBMISSION OPTIONS**

PLEASE NOTE: \*ALL RECEIVED FAXES AND EMAIL WILL BE SENT A CONFIRMATION

① FAX - 1.855.381.2157


② EMAIL - customerservice@parisorthotics.com

**DISCLAIMER:**

Colour of materials may change without notification.

**SHIP VIA**

- Include this order with my next orthotic shipment at no charge
- Ship ASAP via UPS - bulk shipping charge will be added to your order

**SIZING & UNITS**

**SIZING**
**SOFT**
**GENTLE**

GUIDE	SIZING		QUANTITY	PRICE PER PAIR	QUANTITY	PRICE PER PAIR
	MENS	WOMENS				
XS	4 - 5	5 - 7		PARIS ORTHOTICS LAB SERVICES ONLY		RETAIL CONSUMERS CONTACT YOUR CLINICIAN
SM	6 - 7	8 - 9				
M	8 - 9	10 - 11				
LG	10 - 11	12 - 13				
XL	12 - 13	14 +				
XXL	14 - 15	-				
<b>SUBTOTALS</b>			# _____	\$ _____	# _____	\$ _____

**TOTAL** (Before Taxes & Shipping)

**VOLUME DISCOUNT** - Orders over \$200 are eligible for 10% discount