



CARBON FIBRE SPRING PLATE

INTERNAL SHOE STIFFENER

A simple and effective solution for limiting motion across the forefoot and aiding in propulsion

FEATURES

- Extremely thin, light and durable made from premium quality pre-preg carbon
- Pre-shaped left or right and available in 4 sizes
- Contoured design to provide ideal toe spring and heel elevation for most athletic, walking and orthopaedic footwear designs
- Excellent interface with custom foot orthotics
- A practical alternative to some common shoe modifications

CLINICAL INDICATIONS

- Post-surgical care option when transitioning from post-op shoe to regular footwear after bunion, hammertoe and other forefoot surgeries
- Management solution for turf toe, forefoot stress fractures, sesamoiditis and Lisfranc injuries
- Protection and toe off assistance for pathologies affecting the midfoot and forefoot including rheumatoid arthritis, hallux rigidus or limitus, Freiberg's disease, diabetes and forefoot amputations



Sizing: Plates are available left and right in 4 sizes



ORTHOTICS

CARBON FIBRE SPRING PLATE

SHIP TO

ACCOUNT NAME _____

ACCOUNT # _____

CLINICIAN _____

ADDRESS _____

CITY _____

PROVINCE _____ POSTAL CODE _____

PHONE _____ FAX _____

EMAIL _____

CONSIDERATIONS

- Maximum body weight tolerance is 280 lbs.
- Important to use in shoes with adequate forefoot rocker profile.
- Plates are not returnable for credit however size exchange requests will be accommodated for **unused and unmodified plates only**.

PRODUCT GUARANTEE

The Paris Carbon Fibre Spring Plate is guaranteed against breakage for 6 months. Grinding, heating, and reshaping risks delamination/fracture and will nullify warranty.

ORDER SUBMISSION OPTIONS

PLEASE NOTE: *ALL RECEIVED FAXES AND EMAIL WILL BE SENT A CONFIRMATION

① FAX - 1.855.381.2157

② EMAIL - customerservice@parisorthotics.com

BILL TO

ORDERS PLACED WITHOUT AN ACCOUNT # ARE PAYABLE BY CREDIT CARD.
PROVIDE PAYMENT INFORMATION BELOW OR CALL CUSTOMER SERVICE TO PAY-BY-PHONE.

CARD NO. _____ - _____ - _____ - _____

CARD TYPE VISA ☐ MASTERCARD ☐ AMEX ☐

EXPIRY _____ / _____

SHIP VIA

- ☐ Include this order with my next orthotic shipment at no charge
- ☐ Ship ASAP via UPS - bulk shipping charge will be added to your order

ORDER DETAILS - SIZING & UNITS

*Orders of 8 or more units are eligible for a \$5.00 discount per unit.
Recommended **STARTER KIT**: one of each size for Left & Right.

SIZE	DEVICE LENGTH tolerance +/- 2mm		PRICE / UNIT	*BULK ORDER PRICE / UNIT	LEFT / RIGHT	QUANTITY	ITEM TOTAL (\$)
XS	239 MM	INSOLE MEASUREMENT			L		
S	259 MM				L		
M	279 MM				L		
L	298 MM				L		
XS	239 MM				R		
S	259 MM				R		
M	279 MM				R		
L	298 MM				R		
Before taxes and shipping							