

Paris

ORTHOTICS

# CARBON FIBRE SPRING PLATE

INTERNAL SHOE STIFFENER

A simple and effective solution for limiting motion across the forefoot and aiding in propulsion

## FEATURES

- Extremely thin, light and durable made from premium quality pre-preg carbon
- Pre-shaped left or right and available in 4 sizes
- Contoured design to provide ideal toe spring and heel elevation for most athletic, walking and orthopaedic footwear designs.
- Excellent interface with custom foot orthotics
- A practical alternative to some common shoe modifications

## CLINICAL INDICATIONS

- Post-surgical care option when transitioning from post-op shoe to regular footwear after bunion, hammertoe and other forefoot surgeries
- Management solution for turf toe, forefoot stress fractures, sesamoiditis and Lisfranc injuries
- Protection and toe off assistance for pathologies affecting the midfoot and forefoot including rheumatoid arthritis, hallux rigidus or limitus, Freiberg's disease, diabetes and forefoot amputations



**Sizing:** Plates are available left and right in 4 sizes



ORTHOTICS

# CARBON FIBRE SPRING PLATE

## SHIP TO

ACCOUNT NAME \_\_\_\_\_

ACCOUNT # \_\_\_\_\_

CLINICIAN \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY \_\_\_\_\_

PROVINCE \_\_\_\_\_ POSTAL CODE \_\_\_\_\_

PHONE \_\_\_\_\_ FAX \_\_\_\_\_

EMAIL \_\_\_\_\_

## CONSIDERATIONS

- Maximum body weight tolerance is 280 lbs.
- Important to use in shoes with adequate forefoot rocker profile.
- Plates are not returnable for credit however size exchange requests will be accommodated for **unused and unmodified plates only**.

## PRODUCT GUARANTEE

This Paris Carbon Fibre Spring Plate is guaranteed against breakage for 6 months.

- Minor grinding (forefoot shape and/or length) will not affect durability or nullify warranty.
- Heating and reshaping will cause delamination and will nullify warranty.

## BILL TO

ORDERS PLACED WITHOUT AN ACCOUNT # ARE PAYABLE BY CREDIT CARD. PROVIDE PAYMENT INFORMATION BELOW OR CALL CUSTOMER SERVICE TO PAY-BY-PHONE.

CARD NO. \_\_\_\_\_ - \_\_\_\_\_ - \_\_\_\_\_ - \_\_\_\_\_

CARD TYPE VISA  MASTERCARD  AMEX

EXPIRY \_\_\_\_\_ / \_\_\_\_\_

## ORDER SUBMISSION OPTIONS

PLEASE NOTE: \*ALL RECEIVED FAXES AND EMAIL WILL BE SENT A CONFIRMATION

- ① FAX - 1.855.381.2157
- ② EMAIL - customerservice@parisorthotics.com

## SHIP VIA

- Include this order with my next orthotic shipment at no charge
- Ship ASAP via UPS - bulk shipping charge will be added to your order

## ORDER DETAILS - SIZING & UNITS

\*Orders of 8 or more units are eligible for a \$5.00 discount per unit. Recommended **STARTER KIT**: one of each size for Left & Right.

SIZE	DEVICE LENGTH tolerance +/- 2mm	PRICE / UNIT	*BULK ORDER PRICE / UNIT	LEFT / RIGHT	QUANTITY	ITEM TOTAL (\$)
XS	239 MM	<b>INSOLE MEASUREMENT</b>		L		
S	259 MM			L		
M	279 MM			L		
L	298 MM			L		
XS	239 MM			R		
S	259 MM			R		
M	279 MM			R		
L	298 MM			R		

Before taxes and shipping