



# THE STRASSBURG SOCK

## USER-FRIENDLY NIGHT SPLINT

**THE STRASSBURG SOCK IS A VERY EFFECTIVE, SURPRISINGLY COMFORTABLE AND COST EFFECTIVE COMPLIMENT TO TRADITIONAL PLANTAR FASCIITIS TREATMENT REGIMES.**

Conventional hard plastic night splints maintain tension on the fascia by holding the ankle at or near a neutral position. The Strassburg sock provides a more complete and specific tension on the plantar fascia by passively dorsiflexing the toes (and engaging the windlass mechanism) as well as by holding the ankle in neutral position. Used in combination with custom foot orthotics, anti-inflammatory and other conservative therapies, the Strassburg sock night splint significantly improves clinical outcomes.

### ADVANTAGES

- A breathable and washable tubular knit “sock” design that is comfortable enough to sleep in
- A unique hook & lock adjustable strap design that dorsiflexes the toes as well as holds the foot in a neutral position
- No restriction of the ankle joint

### SIZE OPTIONS

- Adult Regular ( up to 16” calf diameter)
- Adult Large ( > 16” calf diameter)
- A pediatric version is also available for use in treating ankle equinus



## THE STRASSBURG SOCK

### SHIP TO

ACCOUNT NAME \_\_\_\_\_

ACCOUNT # \_\_\_\_\_

CLINICIAN \_\_\_\_\_

ADDRESS \_\_\_\_\_

CITY \_\_\_\_\_

PROVINCE \_\_\_\_\_ POSTAL CODE \_\_\_\_\_

PHONE \_\_\_\_\_ FAX \_\_\_\_\_

EMAIL \_\_\_\_\_



### ORDER SUBMISSION OPTIONS

PLEASE NOTE: \*ALL RECEIVED FAXES AND EMAIL WILL BE SENT A CONFIRMATION

① FAX - 1.855.381.2157

② EMAIL - customerservice@parisorthotics.com

### SHIP VIA

- Include this order with my next orthotic shipment at no charge
- Ship ASAP via UPS - bulk shipping charge will be added to your order

### \* SPECIAL ORDER ITEM

 The Pediatric Sock is a **SPECIAL ORDER** item.  
 Shipping turnaround times average 10-14 business days

## SIZING & UNITS

QUANTITY	SIZE	\$ PER UNIT	ITEM TOTAL
	Regular (Up to 16" Calf Diameter )	PARIS ORTHOTICS LAB SERVICE CUSTOMERS ONLY	
	Large (> 16" Calf Diameter )		
	* Paediatric (Up to 12" Calf Diameter )		

### VOLUME DISCOUNT - Price Breaks

6 - 12	(assorted variety).....	5%
13 - 18	(assorted variety).....	10%
19+	(assorted variety).....	15%

**DISCOUNT**
**SUB TOTAL**
*plus applicable taxes*