

\* REQUIRED FIELD

DATE

MM / DD / YYYY

\* ACCOUNT NAME/No. \_\_\_\_\_

\* PATIENT NAME FIRST \_\_\_\_\_ LAST \_\_\_\_\_

\* SHOE SIZE \_\_\_\_\_

\* CLINICIAN NAME FIRST \_\_\_\_\_ LAST \_\_\_\_\_

D.O.B. MM / DD / YYYY \_\_\_\_\_

\* WEIGHT \_\_\_\_\_  LB  KG

US-M  UK

PHONE \_\_\_\_\_ FAX \_\_\_\_\_

FOOTWEAR TYPE \_\_\_\_\_

US-W  EU

# 1 \* CHOOSE A DEVICE

## FUNCTIONAL DEVICES

- STANDARD
- DIRECT MILLED (polypropylene only)

## DRESS DEVICES

- STANDARD SLIMLINE  COBRA (polypro only)
- LOW HEEL CUP  CASUAL
- FLAT HEEL CUP \_\_\_\_\_ CM HEELED SHOES

## SPORT DEVICES

- STANDARD
  - Neoprene to toes  Vinyl to mets
- LOW PROFILE
- DIRECT MILLED (polypropylene only)
  - Neoprene to toes  Vinyl to mets
- IMPACT
- SKI  Alpine  Nordic  Snowboard
- SKATE  Hockey  Figure

## MOLDS

- STANDARD
- LOW PROFILE

## SPECIALIZED DEVICES

- Diabetic  Rheumatoid Arthritic
- RA FLEXIBLE MOLD**  
1.5mm Puff, 1.5mm Poron, Rubberflex & Thermocork, Poron Fill, 1.5mm Nyplex
- D STANDARD TRIDENSITY**  
3.0mm P-Cell, 3.0mm Poron, Rubberflex & 50 durometer Fisher Foam, 1.5mm Cushion Cork
- ADULT ACQUIRED FLATFOOT (inclusive of Posterior Tibial Tendon Dysfunction)
- EVA
- UCBL
- ROBERTS WHITMAN
- GAIT PLATE
  - Induce out-toeing  Induce in-toeing

## SUBSEQUENT PAIR ORDER

- FROM WO# \_\_\_\_\_
- DUPLICATE device (no change)
  - CHANGED device (select device and options)
- NUMBER OF PAIRS \_\_\_\_\_
- 3 DAY RUSH (\$20.00)

# 2 OPTIONAL SPECIFICATIONS

## RIGIDITY

- Rigid
- Semi-Rigid
- Semi-Flexible

## SHELL MATERIALS

- Polypropylene
- Performance Rx
- XT Sprint
- TL2100
- Copolymer

## FILL MATERIALS

- Poron
- EVA
- Thermocork

## FOREFOOT EXTENSIONS

### MATERIAL

- Poron  1.5mm  3.0mm
- Nyplex  1.5mm  3.0mm
- Cushion cork  1.5mm  3.0mm
- Puff  1.5mm  3.0mm

### LENGTH

- Sulcus  Toes

## TOP COVER SELECTIONS

### MATERIAL

- Vinyl
- Neoprene 3.0mm
- Micro-suede
- Leather
- Puff  1.5mm  3.0mm
- P-Cell  1.5mm  3.0mm
- \_\_\_\_\_

### COMBINE MATERIAL SELECTED ABOVE WITH:

- Poron  1.5mm  3.0mm

### LENGTH

- Mets  Sulcus  Toes

## BOTTOM COVER SELECTIONS

### MATERIAL

- Agoflex
- Vinyl
- Puff 1.5mm
- Nyplex 1.5mm
- Micro-suede

# 3 POSTING, MODIFICATIONS & ADDITIONS

## POSTING

- Post to calcaneal vertical  Neutral (as they sit)

REARFOOT Post Type  intrinsic  extrinsic  EVA  
 dual density EVA  
 L \_\_\_\_\_°  VR  VLG R \_\_\_\_\_°  VR  VLG

FOREFOOT Post Type  intrinsic  extrinsic L  R  tip L  R  sulcus  
 L  R  corner L  R  toes  
 L \_\_\_\_\_°  VR  VLG R \_\_\_\_\_°  VR  VLG

Skive L \_\_\_\_\_ mm skive R \_\_\_\_\_ mm skive

Inversion L \_\_\_\_\_° inverted R \_\_\_\_\_° inverted

Motion L \_\_\_\_\_° of motion R \_\_\_\_\_° of motion

## CAST DRESSING

- Minimum  Moderate  Maximum  \_\_\_\_\_

## SHELL MODIFICATIONS

- L  R  Heel Cup Depth  12mm  14mm  16mm  18mm  20mm
- L  R  1st Met Cut Out  with extrinsic tip (support) post
- L  R  1st Ray Cut Out
- L  R  Fascial Accommodation (not available in Direct Milled)
- L  R  Medial Flange
- L  R  Lateral Flange
- L  R  Navicular/Cuneiform Sweet Spot w/ medial flange - AS MARKED
- L  R  Heel Spur Accommodation (not available in Direct Milled)
- L  R  Heel Hole  with poron plug
- L  R  Rigid 1st Extension  sulcus  toes
- L  R  Orthotic Width  wide  narrow

## ADDITIONS

- L  R  Heel Lift  EVA  poron  loose \_\_\_\_\_ mm
- L  R  Heel Pad  1.5mm  3.0mm
- L  R  Horseshoe Spur
- L  R  Metatarsal Pad
- L  R  Metatarsal Bar
- L  R  Metatarsal Raise
- L  R  Neuroma Pad
- L  R  Morton's Extension
- L  R  Reverse Morton's Extension  cushion cork  poron  EVA
- L  R  Functional Hallux Limitus Accommodation
- Lesion Accommodation L  1  2  3  4  5 R  1  2  3  4  5
- L  R  Amputation Accommodation  partial  full

# 4 SPECIAL INSTRUCTIONS

\_\_\_\_\_

\_\_\_\_\_

\_\_\_\_\_



